

Heart rate:	69
P duration:	- ms
PR:	- ms
QRS dur.:	161 ms
T peak:	332 ms
QT:	460 ms
QTc:	495 ms
RR Interval:	864 ms



normal PathECG	
QRS:	23 %
ST:	0 %
T-wave:	0 %

normal WaveECG	
QRS:	27 %
ST:	0 %
T-wave:	0 %

Look on next pages for more details


Patient ID: CienECG examples - ECG : LBBB female 82 id-15906.ecg

Report number: 0001000100000131



normal WaveECG

QRS:	27 %
ST:	0 %
T-wave:	0 %





normal PathECG

QRS:	23 %
ST:	0 %
T-wave:	0 %

Patient ID: CienECG examples - ECG : LBBB female 82 id-15906.ecg

Report number: 0001000100000131



This version of CineECG is a research only beta version. CineECG is a brand of ECG Excellence BV, The Netherlands. CineECG analysis is based on patented technology (Patents pending). CineECG does not provide a diagnostic analysis. CineECG shows: a) the deviation between normal expected CineECG patterns and the CineECG derived from 12 lead ECG data; b) the deviation between normal ECG amplitudes and the median beat ECG. CineECG uses electro-anatomical modelling and proprietary algorithms. CineECG results can be influenced by the quality of the patient recorded ECG data (influenced by ECG electrodes positions on the torso, the orientation and rotation of the heart, the quality of the ECG signal recording device and other characteristics outside the control of ECG Excellence).